

Evaluation of 4 Wheeled Walkers

		PROMEDICS PR30283/30	PROMEDICS PR30283/36	Days 106	Days 105
	Castor	204mm	165mm	200mm	200mm
	Product Weight	7.6kg	11kg	8.3kg	7kg
	Seat Height	584mm	584mm	625	540
	Seat Width	300mm	30"	360	360
	Seat Depth	300mm	29"	330	330
	Approx Cost	£48.32	£48.32	£42.94	£35.47
1.	Height adjustability	841 951mm	830 – 940 good	870 990mm	825 – 945 Good range of heights
2.	Ease of adjustment	Screw	Ok		Easily altered by therapist
3.	Total width open Depth	620 – 710mm	23"	610 – 640mm	590mm
4.	Internal width	440cm	19"	46cm	450mm
5.	Ease of folding e.g. catch	Basket has to be removed. Thin strap.	Poor – replace shoelace	Bag good	Folds easily
6.	Brake action	Poor on strap	Poor moves on sit/stand	Good	Easy to put on/off
7.	Weight limit	114kg	114kg	125kg	125kg 18 stone
8.	Sturdiness		Poor – safety brace not locking after folding	√	Not easy to tip
9.	Manoeuvrability	Ok	Good	√	Performed well on a slope
10	Positive comments	Stick/umbrella holder	Good manoeuvrability	Higher version for taller people	Not easy to tip Bag not basket Cover over backrest
11	Negative comments	Soft foam back. Basket	Backrest due to being shaped, not suitable for everyone		Backrest foam (under cover). May deteriorate quickly.
12	Any other comments		Instructions for use not clear after 6		Brakes seemed quite good. Has a rubber section between metal and metal
13	Overall preference				Preferred Model Add to contract

Evaluation of 4 Wheeled Walkers

		Roma 2463/r	Roma 2462/r
	Castor	200mm	150mm
	Product Weight	7kg	6.5kg
	Seat Height	550mm	550mm
	Seat Width	310mm	310mm
	Seat Depth	320mm	320mm
	Approx Cost	£39.27	£34.94
1.	Height adjustability	810 - 920mm good	810 – 920mm
2.	Ease of adjustment	Ok	
3.	Total width open Depth	65cm	600mm
4.	Internal width	19.5”	42 – 69
5.	Ease of folding e.g. catch	Good	Good
6.	Brake action	Moving when brakes on all surfaces	Good
7.	Weight limit	114kg	114kg
8.	Sturdiness	Not good due to brakes	Ok
9.	Manoeuvrability	Ok	Good
10.	Positive comments	Has bag not basket so easier to fold	Has bag not basket so easier to fold
11.	Negative comments	Backrest not durable as soft foam	Soft back
12.	Any other comments		
13.	Overall preference		

Days 105 agreed by Mobility Sub Group as the preferred contract product