

Suitable for Corner Baths!


The AquaJoy Premier Plus Bathlift

The AquaJoy Premier Plus and Junior Plus Bathlifts are suitable for all types of Baths including Corner Baths without any modification to the bathlift!

The simplicity of the AquaJoy Bathlift means there is no need to modify the bathlift in any way.

Simply position the bathlift so that the side flap is over the lip of the bath on the transferring side for safe transportation on and off the bathlift. If you have dimples or bobbles on your bath then extra middle suction feet should be added.

Clients using the AquaJoy Premier Plus Bathlift in a corner bath will feel safe due to the Strength and Stability of the product.
(26.7 stone weight limit)



AquaJoy Bathlifts Ltd . Unit 52 . Consett Business Park
Consett . Co. Durham . DH8 6BN
Tel: +44(0)1207 501 555 . Fax: +44(0)1207 599 789
E-mail: info@aqua-joy.com . Web: www.aqua-joy.com

Contact Rebecca
+44(0) 1207 501 555
rgraham@aqua-joy.com