

Dali Bed Range - Supplier Lisclare

Representative: Colin Frankland

Tele: 08708502384

If a cot is ever being considered, it is important to bear in mind that the Dali Lisclare contract bed base can be used and all extra parts with which to make the bed into a cot can be added. This means that as the child grows it can be turned into a bed and extended if need be and greatly reduces the price of the overall cot.

General Background Information

Lisclare dali beds are a family of beds where all parts are common and interchangeable.

There is a low entry version available 8" off the ground also there is a 180cm short version for children.

An extra extension piece can be added to this, to extend the cot into a bed, allowing for growth.

A point worth noting - cots do take up a lot of valuable space within a warehouse, the Lisclare cot parts can be stored flat and built up as needed.

Oasis Cot

Cost – with padding approx £3k without padding £2,300

- Door entry system allows instant access
- Previous model had keys at the top of the doors, which were situated on the outside, to open them. They now have a slide mechanism at the top, which shuts the doors securely. There are extra locks at the base of the doors, which can be operated via a lever handle. The doors can be open partly or swung completely out.
- The wooden bars can be taken off and the sides replaced with metal bed rails if required
- All handsets have lock out facilities
- The wire lead used for the handset is reinforced
- If children are able to pull the velcro padding off, then bolts can be added to secure the padding
- The sides can be made to whatever height that is requested
- A hole can be made in the side of the frame or padding for a feeding tube if required

- The size of the mattress needs to be considered before placing an order. If padded sides are required within the cot, then the mattress dimensions need to be smaller
- Smaller beds can be made for confined spaces
- The lead time for product manufacture is 3 – 4 weeks
- All products made in Wrexham, UK

Bariatric Bed

The company make 3 bariatric beds in 90cm, 100cm or 120cm width. The 100cm width is proving to be the most popular.

Cost - approx £1245, with mattress an extra £1000

Beds

The company make a 200cm length adult size dali bed, which can be extended to 220cm with an extension piece. This bed is on the Derbyshire contract.

There is a 180cm child sized bed which is the base for the Oasis Cot and an extension can be added as the child grows making it a 200cm length.

NB only one extension piece can be added to a bed frame

Sometimes a wider bed may be needed even if it is for a child, for example a child may be very strong and may require a wider bed to accommodate being physically active.

A padded surround can be used without wooden sides on the dali bed

7.1.10

Judy Green
Lead Practitioner
SROT
Derbyshire ICES

17/12/2009